

# Polarimeter Updates

- Scaler readouts
- Running DAQ
- Pedestal adjustment
- Missing channel
- DATA!!

April 29, 2015  
Kei Moriya, Michael Dugger

# Translation Table

- Convert from rocid, slot, channel to detector
- Enables grabbing objects with JANA framework
- Also necessary for scalers to be read out
- Within ccdb, set run > 2881, is Translation/DAQ2detector
- Run ccdb -r 2882 dump Translation/DAQ2detector to see contents
- Thanks to Paul M. for putting these in

# Translation Table

- Script used is in `online/packages/TranslationTable/tt_add_tpol_sectors.py`
- Check what is in `tt.db`:

```
> sqlite3 tt.db
```

```
> sqlite> .headers on
```

```
sqlite> .schema
```

```
CREATE TABLE BCAL(module INT, layer INT, sector INT, end TEXT, adc_chanid INTEGER, tdc_chanid INTEGER, disc_chanid INTEGER);
CREATE TABLE CDC(ring INT, straw INT, adc_chanid INTEGER);
CREATE TABLE Channel(chanid INTEGER PRIMARY KEY, moduleid INTEGER, name TEXT, channel INT, system TEXT, col_name TEXT, enable INTEGER DEFAULT 1);
CREATE TABLE Crate(crateid INTEGER PRIMARY KEY, name TEXT, area TEXT, rack INT, location TEXT, SN TEXT, Function TEXT, rocid INT, host TEXT, IP TEXT);
CREATE TABLE Detector_Hierarchy( id INTEGER PRIMARY KEY, parent_id INTEGER, name TEXT, type TEXT, chanid INTEGER, mtime DATETIME DEFAULT CURRENT_TIMESTAMP);
CREATE TABLE FCAL(row INT, col INT, adc_chanid INTEGER);
CREATE TABLE FDC_Cathodes(package INT, chamber INT, view TEXT, strip INT, strip_type TEXT, adc_chanid INTEGER);
CREATE TABLE FDC_Wires(package INT, chamber INT, wire INT, tdc_chanid INTEGER, disc_chanid INTEGER);
CREATE TABLE Module(moduleid INTEGER PRIMARY KEY, crateid INTEGER, slot INT, type TEXT, SN TEXT);
CREATE TABLE PS(side TEXT, id INT, adc_chanid INTEGER);
CREATE TABLE PSC(id INT, adc_chanid INTEGER, tdc_chanid INTEGER, disc_chanid INTEGER);
CREATE TABLE RF(system TEXT, type TEXT, tdc_chanid INT, adc_chanid INT);
CREATE TABLE ST(sector INT, adc_chanid INTEGER, tdc_chanid INTEGER, disc_chanid INTEGER);
CREATE TABLE TAGH(id INT, adc_chanid INTEGER, tdc_chanid INTEGER, disc_chanid INTEGER);
CREATE TABLE TAGM(col INT, row INT, adc_chanid INTEGER, tdc_chanid INTEGER, disc_chanid INTEGER);
CREATE TABLE TOF(plane INT, bar INT, end TEXT, adc_chanid INTEGER, tdc_chanid INTEGER, disc_chanid INTEGER);
CREATE TABLE TPOL(sector INT, ring INT, adc_chanid INT);
sqlite> .exit
```

# Translation Table TPOL

```
<crate number="84" type="VXS">
  <slot number="1" type="CPU">
  </slot>
  <slot number="11" type="CTP">
  </slot>
  <slot number="12" type="SD">
  </slot>
  <slot number="13" type="FADC250">
    <channel number="0" detector="TPOL" sector="1" ring="0"/>
    <channel number="1" detector="TPOL" sector="2" ring="0"/>
    <channel number="2" detector="TPOL" sector="3" ring="0"/>
    <channel number="3" detector="TPOL" sector="4" ring="0"/>
    <channel number="4" detector="TPOL" sector="5" ring="0"/>
    <channel number="5" detector="TPOL" sector="6" ring="0"/>
    <channel number="6" detector="TPOL" sector="7" ring="0"/>
    <channel number="7" detector="TPOL" sector="8" ring="0"/>
    <channel number="8" detector="TPOL" sector="9" ring="0"/>
    <channel number="9" detector="TPOL" sector="10" ring="0"/>
    <channel number="10" detector="TPOL" sector="11" ring="0"/>
    <channel number="11" detector="TPOL" sector="12" ring="0"/>
    <channel number="12" detector="TPOL" sector="13" ring="0"/>
    <channel number="13" detector="TPOL" sector="14" ring="0"/>
    <channel number="14" detector="TPOL" sector="15" ring="0"/>
    <channel number="15" detector="TPOL" sector="16" ring="0"/>
    <channel number="16" detector="RF" system="PS" type="FADC250"/>
  </slot>
```

```
  <slot number="14" type="FADC250">
    <channel number="0" detector="TPOL" sector="17" ring="0"/>
    <channel number="1" detector="TPOL" sector="18" ring="0"/>
    <channel number="2" detector="TPOL" sector="19" ring="0"/>
    <channel number="3" detector="TPOL" sector="20" ring="0"/>
    <channel number="4" detector="TPOL" sector="21" ring="0"/>
    <channel number="5" detector="TPOL" sector="22" ring="0"/>
    <channel number="6" detector="TPOL" sector="23" ring="0"/>
    <channel number="7" detector="TPOL" sector="24" ring="0"/>
    <channel number="8" detector="TPOL" sector="25" ring="0"/>
    <channel number="9" detector="TPOL" sector="26" ring="0"/>
    <channel number="10" detector="TPOL" sector="27" ring="0"/>
    <channel number="11" detector="TPOL" sector="28" ring="0"/>
    <channel number="12" detector="TPOL" sector="29" ring="0"/>
    <channel number="13" detector="TPOL" sector="30" ring="0"/>
    <channel number="14" detector="TPOL" sector="31" ring="0"/>
    <channel number="15" detector="TPOL" sector="32" ring="0"/>
  </slot>
  <slot number="21" type="TI">
  </slot>
</crate>
```

RF for PS was there already, need to remove later

# Scalars

- Monitor rates online
- In future will be monitored by GUI (Paul M.) and saved in archiver
- For now, log on to rocps2 from hallgw
- Run `> /home/pmatt/controls/epics/R3-14-12-3-RHEL5/drivers/bin/linux-x86/testBoards 1 1`
- 1st arg is to use raw values (0) or rates (1),  
2nd arg is debug flag
- Slots 13 and 14 are triplet sectors
- Need to add in rings later

```
FADC250: Slot, channel, time, scalars = 13, 1, 0.97856, 0
FADC250: Slot, channel, time, scalars = 13, 2, 0.97856, 0
FADC250: Slot, channel, time, scalars = 13, 3, 0.97856, 0
FADC250: Slot, channel, time, scalars = 13, 4, 0.97856, 0
FADC250: Slot, channel, time, scalars = 13, 5, 0.97856, 0
FADC250: Slot, channel, time, scalars = 13, 6, 0.97856, 0
FADC250: Slot, channel, time, scalars = 13, 7, 0.97856, 0
FADC250: Slot, channel, time, scalars = 13, 8, 0.97856, 0
FADC250: Slot, channel, time, scalars = 13, 9, 0.97856, 0
FADC250: Slot, channel, time, scalars = 13, 10, 0.97856, 0
FADC250: Slot, channel, time, scalars = 13, 11, 0.97856, 0
```

# Running DAQ

- Sasha has set up fADCs and DAQ so we can take self-triggered data
- Log in to hdops, run `rcm.sh`; this opens the DAQ GUI
- Configuration files are in `/gluonfs1/home/hdops/CDAQ/config/hd_psc/`
- Files `tpol_pulser.conf` and `tpol_self.conf`
- Final part is for TPOL settings, still need to be set to appropriate values





# Pedestal Adjustment

- Mark Dalton showed how to adjust pedestals for each channel
- Within config file set DAC values which move the pedestals
- Standard is to have pedestals near 100 ADC counts (max count 4095)

note log scale

For each channel pedestal is  
cycled through values of 0 - 1000  
and scaler rates are read out





channels 0 - 15







# Channels 0 - 15

- Pedestals now centered at ADC values of 100

before



after







→ ADC counts







# Channels 16 - 31

- Pedestals now centered at ADC values of 100
- Ch 31 had very large (~50 mV) baseline fluctuation on scope

before



after





→ ADC counts

# Pedestal Adjustment

- Fit pedestal distributions with Gaussian functions


before adjustment



Width is much wider compared to other detectors,  
probably due to baseline 60Hz noise + fluctuation

# Pedestal Adjustment



- Fit pedestal distributions with Gaussian functions



- Ch 15 was not adjusted yet (next slides)
- Ch 16 and 19 may not be working (narrowness of distribution suggests this)
- Ch 31 has baseline fluctuations

# Missing Channel

- Pointed out by Sasha that one channel was missing (slot 13, channel 15)
- Confusion in electronics group about where to put RF cable
- Went in during controlled access, plugged cable





# Other Checks

- Checked what I could
- HV, LV seem fine, fans are working
- Not sure which vacuum level is for polarimeter




# DATA!!!!

- Sasha started configuring the settings and taking data on Monday night
- Initial look is run 2941, ~20k events
- Grabbing fADC250 pulse integral objects without JANA framework, will work on this when I have time
- Data is from PS trigger (coincidence in left and right arms of coarse PS), 100 samples are read out

# DATA!!!!





- Raw waveforms from hit (100 samples, 400 ns)
- Decay time is  $\sim 20$  ns, much shorter than what we are used to

Run 2914



# DATA!!!!

Run 2914





# More Data

1. [https://halldweb.jlab.org/cgi-bin/data\\_monitoring/run\\_conditions.pl](https://halldweb.jlab.org/cgi-bin/data_monitoring/run_conditions.pl)

- PS trigger data from night of April 28 - 29


run	# events <sup>1</sup>	TPOL converter	radiator
2914	31.8k	none	
3007	202.8k	none	
3012	201.1k	none	
3015	30k	25 $\mu$ m	
3016	8k	50 $\mu$ m	
3017	7.5k	50 $\mu$ m	
3018	2.69M	none	diamond
3019	1.81M	none	diamond
3020	2.18M	none	amorphous
3021	2.34M	none	amorphous
3022	3.42M	none	amorphous
3023	183k	none	amorphous
3024	107k	none	amorphous

# Analysis

- No qualitative difference between runs with and without converter
- Will work on getting TPOL into GlueX framework so that we can relate hits with other systems (PS, tagger)
- Need to understand features


# # of Hits/Event

- Run 3015
- Same sector can have multiple hits
- Sasha may be trying to read out all channels regardless of pulse heights



# # of Hits By Sector

- Run 3015
- Missing 2 sectors
- Would partly explain tendency to get 30 hits/event if reading out all channels that have ADC counts  $> 110$



# Looking Forward

- Need to check a few things in hall tomorrow
- Analyze current data, more coming - have requested different converter combinations